English	Instructions for use: Dental Implants	3
Español	Instrucciones de uso: Implantes Dentales	6
русский	Инструкция по применению: Зубные имплантаты	10
Deutsch	Gebrauchsanweisung: Dentalimplantate	14
Français	Instructions d'utilisation: Implants Dentaires	18
Italiano	Istruzioni per l'uso: Impianti Dentali	22
Português	Instruções de utilização: Implantes Dentários	25
Türkçe	Kullanım talimatları: Dental İmplantlar	28
简体中文	使用说明: 牙种植体	31
日本語	使用説明書: 歯科インプラント	34
한국어	적응증: 치과 시술용 임플란트	39
عربي	إرشادات الاستخدام: غرسات الأسنان	44
Polski	Instrukcja obsługi: Implanty dentystyczne	47

manufactured by:  
**BioHorizons**  
2300 Riverchase Center  
Birmingham AL, 35244 USA  
TOLL FREE 888.246.8338  
TEL 205.967.7880 FAX 205.870.0304  
www.biohorizons.com

EU Authorised Representative:  
**QUALITY FIRST INTERNATIONAL**  
Suites 317/318 Burford Business Centre  
11 Burford Road, Stratford  
London E15 2ST United Kingdom  
Telephone +44-208-221-2361  
Telefax +44-208-221-1912

Symbol	Symbol Description
	Caution
	Consult Instructions for use
	Manufacturer
	BioHorizons products carry the CE mark and fulfill the requirements of the Medical Devices Directive 93/42/EEC
	Reference/ article number
	Lot/ batch number
	Do not re-use
	Use-by-date
	Sterile by gamma irradiation
	Date of manufacture
<b>Rx Only</b>	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician

## ENGLISH

This document supersedes all prior revisions. Original language is English.

This document applies to BioHorizons Dental Implants.

### DESCRIPTION

BioHorizons implants are manufactured from biocompatible titanium alloy and are single use only. Some BioHorizons implants come packaged with a 3inOne abutment pre-mounted. 3inOne abutments are for single patient use only. This abutment may be used as: (1) an implant mount; (2) a closed tray transfer when used with a Ball-top Screw and (3) a temporary or final abutment. "Peel-and-stick" labels supplied on the implant package contain important product information and should be applied to the patient's record in the event future reference is necessary.

### INDICATIONS FOR USE

BioHorizons implants are intended for use in the mandible or maxilla for use as an artificial root structure for single tooth replacement or for fixed bridgework and dental retention. The implants may be restored immediately (1) with a temporary prosthesis that is not in functional occlusion or (2) when splinted together for multiple tooth replacement or when stabilized with an overdenture supported by multiple implants.

Tapered Short implants may be restored using delayed loading, or with a terminal or intermediate abutment for fixed or removable bridgework, and for overdentures.

Tapered Internal 3.0, Tapered Tissue Level 3.0, and Laser-Lok 3.0 Implants may be used as an artificial root structure for single tooth replacement of mandibular central and lateral incisors and maxillary lateral incisors. The implants may be restored immediately (1) with a temporary prosthesis that is not in functional occlusion, (2) when splinted together as an artificial root structure for multiple tooth replacement of mandibular incisors, or (3) for denture stabilization using multiple implants in the anterior mandible and maxilla. The implants may be placed in immediate function when good primary stability has been achieved and with appropriate occlusal loading.

One-piece 3.0 Implants may be used (1) as an artificial root structure for single tooth replacement of mandibular central and lateral incisors and maxillary lateral incisors. The implant may be immediately restored with a temporary prosthesis that is not in functional occlusion; (2) when splinted together as an artificial root structure for multiple tooth replacement of mandibular incisors. The implants may be restored after a period of time or placed in immediate function; or (3) for denture stabilization using multiple implants in the anterior mandible and maxilla. The implants may be restored after a period of time or placed in immediate function.

Overdenture Implants may be used for denture stabilization using multiple implants in the anterior mandible and maxilla. The implants may be restored after a period of time or placed in immediate function.

### CONTRAINDICATIONS

BioHorizons dental implants should not be used in patients who have contraindicating systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity.

Tapered Internal 3.0, Tapered Tissue Level 3.0, and Laser-Lok 3.0 Implants are contraindicated for single tooth replacement of maxillary central incisors.

One-piece 3.0 Implants are contraindicated for single use to replace maxillary central incisors, nor should they be used as abutments for multiple bridgework cases other than as specified in the Indications section of these instructions. They should not be placed in occlusal function for a minimum of three months post-implantation in high-density bone, or approximately six months in lower-density bone.

### **DIRECTIONS FOR USE**

Proper surgical procedures and techniques are the responsibility of the medical professional. Each dentist must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of postgraduate dental implant education and strict adherence to the instructions and procedures for use that accompany BioHorizons implant products. NOTE: A pre-operative 30-second rinse with a 0.12% Chlorhexidine Digluconate solution is recommended. (The influence of 0.12% Chlorhexidine Digluconate Rinses on the Incidence of Infectious Complications and Implant Success. Lambert, Paul M., et al, J Oral Maxillofacial Surgery 55:25-30, 1997, Suppl5).

Tapered Internal Plus Implants - IMPORTANT! When using the TSK2021 (Tapered Internal Surgical Kit) to place the 3.8mm diameter implants, 3.0mm implant-level drivers are also required (TP31DH, TP31DR).

Tapered Internal 3.0, Tapered Tissue Level 3.0, and Laser-Lok 3.0 Implants - IMPORTANT! (1) The implant procedure manual should be reviewed prior to attempting placement. (2) The osteotomy may be prepared with any BioHorizons 2.0mm and 2.5mm drill used in sequence. (3) Torque of 70Ncm or greater may compromise the mechanical integrity of the implant. In dense (D1-D2) bone further preparation of the site may be required to avoid applying excessive torque to the implant connection. (4) To further prepare the site, the crestal bone drill (TP3CBD) and the bone tap (122-900) may be used. (5) The clinician should verify that the proper instruments are available and at hand prior to attempting placement.

### **WARNINGS AND PRECAUTIONS**

Appropriate training in proper implant surgery technique is strongly recommended prior to implant use. Improper technique can result in implant failure and/or loss of supporting bone. Appropriate x-ray films and/or CT scans should be utilized to determine (1) if adequate bone is available and (2) the location of important anatomical landmarks, such as the mandibular canal, maxillary sinuses, and adjacent teeth.

Clinician judgment should be used in determining the minimum post-implantation time before placing the implants in occlusal function. An adequate number of implants should be used to provide support and to distribute the load to the abutments. For short implants, clinicians should closely monitor patients for any of the following conditions: periimplant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinicians choose a short implant, then clinicians should consider a two-stage surgical approach, splinting the short implant to an additional implant, and placement of the widest possible fixture. Allow longer periods for osseointegration avoid immediate loading. BioHorizons dental implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of BioHorizons dental implants in the MR environment is unknown. Scanning a patient who has this device may result in patient injury. Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. A certain percentage of implants may fail to achieve or to maintain osseointegration, as demonstrated by mobility, and should be removed. Potential causes of abutment fracture include, but not limited to: casting beyond a 30-degree angle, casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive

fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics. Small diameter implants with angled abutments are intended for the anterior region of the mouth and not intended for the posterior region of the mouth due to possible failure of the implant.

One-piece 3.0 Implants should NOT be placed in full occlusion until a minimum of THREE months in situ from the date of placement in bone of high integrity, and a minimum of SIX months in situ in bone of compromised integrity.

Overdenture Implants MUST be carefully evaluated for stability before selecting an appropriate level of retention in mating overdenture attachment inserts, o-rings, or grommets. Less retention rather than more retention is recommended for initial loading. Overdenture implants are designed for denture stabilization; a well-fitting denture with soft-tissue support is highly recommended.

### **COMPLICATIONS AND ADVERSE EFFECTS**

The risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant and/or abutment; (3) loosening of the abutment screw and/or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; (9) perforation of the maxillary sinus; (10) perforation of the labial and lingual plates; and (11) bone loss possibly resulting in revision or removal.

### **HANDLING AND STERILIZATION**

This product is supplied sterile and should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove from the package only after the correct size has been determined and the surgical site has been prepared. Always handle the product with powder-free gloves, and avoid contact with hard objects that may damage the surface. Implants are single use only, and re-use should not be attempted. Following this guidance eliminates the risk of cross-patient contamination from secondary use of this device. BioHorizons assumes no responsibility for attempted re-use or re-sterilization.

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