English	Instructions for use: Custom abutments	3
Español	Instrucciones de uso: Soportes personalizados	7
русский	Инструкция по применению: Индивидуальные абатменты	11
Deutsch	Gebrauchsanweisung: Abutments nach Maß	16
Français	Instructions d'utilisation: Piliers prothétiques personnalisés	20
Italiano	Istruzioni per l'uso: Pilastrini personalizzati	24
Português	Instruções de utilização: Pilares personalizados	28
Türkçe	Kullanım talimatları: Özel Abutmentler	32
简体中文	使用说明：自定义基牙	36
日本語	使用説明書：カスタム アバットメント	40
한국어	적응증: 맞춤형 지대주	44
عربي	إرشادات الاستخدام: دعائم مخصصة	49
Polski	Instructions for use: Custom abutments	54

BioHorizons  
2300 Riverchase Center  
Birmingham AL, 35244 USA TOLL  
FREE 888.246.8338  
TEL 205.967.7880 FAX 205.870.0304  
www.biohorizons.com

Quality First International OÜ  
Laki 30  
12915 Tallinn  
Estonia  
Tel +372-610-4196

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution
	Electronic instructions for use
	Manufacturer
	BioHorizons products carry the CE mark and fulfill the requirements of the Medical Devices Directive 93/42/EEC
	Reference/ article number
	Lot/ batch number
	Do not re-use
	Use-by-date
	Sterile by gamma irradiation
	Date of manufacture
<b>Rx Only</b>	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician
	EU Authorised Representative
	Do not use if package is damaged
	Medical Device
	Non-sterile
	Home position

This document supersedes all prior revisions. Original language is English.

## DESCRIPTION

BioHorizons CAD/CAM titanium blanks (Custom Ti and Custom Milled) and CAD/CAM titanium bases (Titanium Base and Hybrid Base) are intended for the restoration of BioHorizons Tapered Internal and Internal dental implants and Zimmer® Dental ScrewVent® and Tapered ScrewVent® dental implants within the specific indications of each implant system. The label on each prosthetic package contains important product information including an indication that the prosthetic is supplied non-sterile.

## INDICATIONS FOR USE

BioHorizons CAD/CAM Abutments are dental abutments placed onto a dental implant to provide support for dental prosthetic restorations. The abutments include: 1) Titanium abutment blanks with a pre-machined implant connection where the upper portion may be custom-milled in accordance with a patient-specific design using CAD/CAM techniques; and 2) Titanium bases with a pre-machined implant connection upon which a CAD/CAM designed superstructure may be fitted to complete a two-piece dental abutment. The abutments include an abutment screw for fixation to the underlying implant. The abutments may be used for single-unit (single-tooth) or multiple-unit (bridges and bars) restorations and are compatible for use with BioHorizons Internal and Tapered Internal implant systems and Zimmer® Dental Screw-Vent® and Tapered Screw-Vent® implants with 3.5mm, 4.5mm and 5.7mm internal hex-connection mating platform diameters.

All digitally designed abutments and/or copings for use with BioHorizons CAD/CAM Abutments are intended to be sent to a BioHorizons-validated milling center for manufacture. BioHorizons abutments designed using CAD/CAM techniques must fulfill the BioHorizons allowable range of design parameters.

## CONTRAINDICATIONS

BioHorizons prosthetics should not be used in patients who have contraindicated systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity.

## DIRECTIONS FOR USE

A custom abutment design file is created by the customer using a 3Shape abutment library. The abutment design file is converted to a stereolithography (.stl) file using 3Shape software. The .stl file is converted to a numerical control (.nc) file using an appropriate software application (e.g. SUM3D). After the custom abutment design is uploaded, the milling machine software executes the necessary commands to mill the final abutment. The digitally designed file is to be sent to a registered and listed BioHorizons contract manufacturer for milling. Visit [www.biohorizons.com](http://www.biohorizons.com) for a list of validated milling centers. BioHorizons CAD/CAM Abutments are compatible with commercially available dental CAD/CAM systems, such as 3Shape.

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of dental implant courses and strict adherence to the Instructions for Use that accompany BioHorizons products. The Single-stage Hybrid Abutment Base is used for single-unit restorations where the underlying implant is placed supra-crestally and a single-stage protocol is desired. Abutment screws must be torqued to 30 Ncm to prevent screw loosening. BioHorizons Custom Milled Abutments, Hybrid Abutment Bases and 3.0mm Custom Ti Abutments contain a LIGHT BLUE abutment screw specifically designed for use in those abutments. Use of the LIGHT BLUE

abutment screws with incompatible abutments may result in abutment and/or abutment screw failure. Restorations utilizing angled screw channels are only compatible with YELLOW Precision Angled Screws. Angled screw channels can have a maximum angulation of 15° off the vertical axis of the abutment.

Design limitation parameters for the Custom Ti and Custom Milled Abutments are included in the following table.

Feature	Custom Ti Platforms				Custom Milled Platforms	
	Ø3.0mm		Ø3.5mm, Ø4.5mm, Ø5.7mm		Ø3.0mm, Ø3.5mm, Ø4.5mm, Ø5.7mm	
	Min	Max	Min	Max	Min	Max
Wall thickness	0.3mm	3.8mm	0.3mm	3.6mm	0.3mm	3.7mm
Gingival margin diameter	3.0mm	9.7mm	3.0mm	9.7mm	3.3mm, 3.8mm, 4.6mm, 5.8mm respectively	9.5mm
Gingival margin height	0mm	12.1mm	0mm	12.1mm	0.3mm, 0.2mm, 0.1mm, 0.1mm respectively	12.3mm, 12.5mm, 12.6mm, 12.5mm respectively
Length of the abutment post	*	12.1mm	*	12.1mm	*	12.3mm, 12.5mm, 12.6mm, 12.5mm respectively
Taper of the abutment	0°	90°	0°	90°	0°	90°
Angulation	0°	30°	0°	30°	0°	30°

\* Minimum allowable abutment post dimensions are 3.0mm for screw retained abutment and 4.0mm for cement retained abutment.

Titanium base abutments - BioHorizons recommends the use of Metoxit Z-CAD® dental ceramic (or equivalent) to fabricate the custom zirconia superstructure and RelyX Unicem (or equivalent) to affix the zirconia superstructure to the titanium base.

Design limitation parameters for the custom zirconia superstructure portion of the Titanium Base and Hybrid Base Abutments are included in the following table.

Feature	Titanium Base Platforms		Hybrid Base Platforms	
			Ø3.0mm†, Ø3.5mm, Ø4.5mm, Ø5.7mm	
	Min	Max	Min	Max
Wall thickness	0.4mm	**	0.4mm	**
Gingival margin diameter	5.1mm	**	3.8mm, 3.8mm, 4.6mm, 5.8mm respectively & 5.1mm, 6.1mm, 7.4mm respectively for Single-stage	**
Length of the abutment post	3.9mm	**	4.0mm	**
Taper	Variable based on anatomy	**	Variable based on anatomy	**
Angulation of the final titanium base	0°	20°	0°	20°

† Not available for Single-stage.

\*\*Feature maximum sizes are based on patient anatomy and clinician discretion.

Secure Laser-Lok abutments onto the Laser-Lok Protective Sleeve to protect Laser-Lok during handling.

**WARNINGS AND PRECAUTIONS**

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any BioHorizons Instructions for Use (IFU). Clinicians are responsible for understanding the appropriate technical use of BioHorizons prosthetic

components. Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at [www.biohorizons.com](http://www.biohorizons.com). Contact BioHorizons Customer Care or your local representative with any questions you have regarding specific IFU.

Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. Potential causes of abutment fracture include, but are not limited to: casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics.

The BioHorizons CAD/CAM Abutments have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of BioHorizons CAD/CAM Abutments in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

If any modifications are made to the implant/abutment interface, the abutment may not properly interface with the implant. The FDA considers the modifier of the implant/abutment interface a medical device company subject to FDA rules and regulations.

Titanium base abutments - Altering the size or shape of the titanium base abutment may affect proper fit with the restoration. Leaving bonding agent or particulate in the Laser-Lok surface may affect the connective tissue attachment.

Prosthetics are single patient use only. To eliminate the risk of cross-patient contamination re-use should not be attempted. BioHorizons assumes no responsibility for attempted re-use or re-sterilization between patients.

### **COMPLICATIONS AND ADVERSE EFFECTS**

The risks and complications with prosthetic components and implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) implant and/or abutment breakage; (3) abutment screw and/or retaining screw loosening; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) implant loosening requiring revision surgery; (9) maxillary sinus perforation; (10) labial or lingual plate perforation; and (11) bone loss possibly resulting in revision or removal.

### **HANDLING AND STERILIZATION**

Always handle the product with powder-free gloves and avoid contact with hard objects that may damage the surface. After abutment preparation and before final placement in the restorative site, CAD/CAM abutments must be cleaned and sterilized. The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a sterile container using a broad-spectrum cleaning or disinfecting agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations.
- 2) Brush the product to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 3) Thoroughly rinse product under running utility tap water.
- 4) Place product in the sterile container filled with the prepared detergent solution and sonicate for two (2) minutes minimum.

- 5) Thoroughly rinse product under running utility tap water.
- 6) Spray or wipe product with 70% IPA.
- 7) Blot product dry with clean lint free cloth.

For sterilization, place final-prepped and cleaned abutment in an approved sterilization bag or wrap and run through one of the following qualified sterilization cycles:

<b>Sterilization Cycles</b>				
<b>Reference:</b>	ANSI/AAMI TIR12:2010	ANSI/AAMI TIR12:2010	ANSI/AAMI TIR12:2010	UK HTM 01-01 Part C:2016
<b>Type:</b>	Gravity Steam	Gravity Steam	Prevacuum Steam	Prevacuum Steam
<b>Exposure Time and Temperature:</b>	30 minutes at 121°C (250°F)	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
<b>Minimum Dry Time:</b>	30 minutes	30 minutes	20 minutes	20 minutes

It is recommended to include a 30-minute cool-down period before removing the product from the sterilization bag or wrap.

Zimmer®, Dental ScrewVent®, and Tapered ScrewVent® are registered trademarks of Zimmer Dental, Inc.

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©Zimmer، و©Dental ScrewVent، و©Tapered ScrewVent هي علامات تجاية خاصة بشركة Zimmer Dental, Inc.
Zimmer®, dentystyczne ScrewVent® i stożkowe ScrewVent® to zastrzeżone znaki towarowe Zimmer Dental, Inc.